

Dry Blood Spot Punching System

- Sample codes can be traced back
- Support for time-delay automatic or manual punching
- Laser positioning indication of blood spot
- Match any collection cards
- Professional camera visual detection
- Chinese & English interface is customized on demand

Parameters

Model	DY-03
Throughput	1~96 well plate OR 1~24 positions of 2ml centrifuge tube
Puncher Diameter	3mm/6mm (3mm is standard, and the rest are optional)
Plate Position	1 plate position fits 96-well plate, PCR plate, microplate, 2ml centrifuge tube (optional accessories)
Endoscopic Camera	5 million pixels
Puncher Accuracy	±0.01mm
Operating Humidity	20~80%RH
Operating Temperature	10~40 °C
Program Management	Built-in program
Interfaces	3 USB interfaces
Dimensions (L*W*H)	420*320*460mm
Power Supply	A.C. 220 V, 50 Hz
Weight	22.6 kg



Ordering Information

Cat. #	Description	Qty.
DBS30-E	DY-03 Dry Blood Spot Punching System	1 Set/Carton
XSB003	Neonatal Blood Collection Card	100 Pcs/Box